



## Female Luer Connector

IMM has different design female luer connectors for different user preferences:

- Small wing, big wing, no wing
- Different tube OD
- Material EtO compatible unless otherwise noted

### Female Luer Connector

P/N	Description	Tube OD	Wing	Color	Material	Note
A0042	Female luer lock ID2.0mm	2.0	Small	Transparent	PVC	
A0016	Female luer lock ID2.4mm	2.4	Small	Transparent	PVC	
A0139-2B	Female luer lock ID2.4mm	2.4	Small	Transparent blue tint	PVC	
A0154	Female luer lock ID2.4mm	2.4	Small	Orange	PVC	
A0345	Female luer lock ID2.4mm	2.4	Small	Transparent	PVC	
A0366	Female luer lock ID2.4mm	2.4	Small	Transparent light brown tint	PVC	
A0512	Female luer lock ID2.4mm	2.4	Small	Transparent lavender tint	PC	ETO/Gamma
A0349	Female luer lock ID2.4mm	2.4	No	Transparent	PVC	
A0359	Female luer lock ID2.4mm	2.4	No	Yellow	PVC	
A0360	Female luer lock ID2.4mm	2.4	No	Grey	PVC	
A0017	Female luer lock ID2.75mm	2.8	Small	Transparent	PVC	
A0600	Female luer lock ID2.8mm	2.8	Small	Transparent lavender tint	PC	ETO/Gamma
A0435	Female luer lock ID2.8mm	2.8	Small	Transparent	PVC	
A0354	Female luer lock ID2.8mm	2.8	No	Yellow	PVC	
A0398	Female luer lock ID3.0mm	3.0	Big	Transparent	PVC	
A0276	Female luer lock ID3.2mm	3.2	Small	Transparent	PVC	
A0294	Female luer lock ID3.2mm	3.2	Small	Transparent	PVC	
A0591	Female luer lock ID3.2mm	3.2	Big	Transparent	PVC	
A0463	Female luer lock ID3.2mm	3.2	Big	Transparent	PC	
A0591	Female luer lock ID3.2mm	3.2	Big	Transparent	PVC	
A0018	Female luer lock ID3.5mm	3.5	Small	Transparent	PVC	

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: [service@immc.com.tw](mailto:service@immc.com.tw)

*All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.*

Jan 2021

Female Luer Connector						
P/N	Description	Tube OD	Wing	Color	Material	Note
A0518	Female luer lock ID3.5mm	3.5	Small	Blue translucent	PVC	
A0519	Female luer lock ID3.5mm	3.5	Small	Orange translucent	PVC	
A0204	Female luer lock ID3.6mm	3.6	Small	Transparent	PVC	
A0510	Female luer lock ID3.6mm	3.6	Big	Transparent	PC	
A0588	Female luer lock ID3.8mm	3.8	Small	Transparent	PVC	
A0019	Female luer lock ID4.1mm	4.1	Small	Transparent	PVC	
A0142-2B	Female luer lock ID4.1mm	4.1	Small	Transparent blue tint	PVC	
A0207	Female luer lock ID4.1mm	4.1	Small	Transparent lavender tint	PC	ETO/Gamma
A0343	Female luer lock ID4.1mm	4.1	Small	Transparent	PVC	
A0428	Female luer lock ID4.1mm	4.1	Small	Transparent lavender tint	PC	ETO/Gamma
A0339	Female luer lock ID4.1mm	4.1	Big	Transparent	PVC	
A0440	Female luer lock ID4.1mm	4.1	Big	Transparent lavender tint	PC	ETO/Gamma
A0358	Female luer lock ID4.8mm	4.8	Small	Transparent	PVC	
A0589	Female luer lock ID4.9mm	4.9	Small	Transparent	PVC	
A0211	Female luer lock ID5.3mm	5.3	Small	Transparent	PVC	
A0094	Double female luer lock	n/a	No	Original	PP	

**Innovative Medical Manufacturing Company**

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

*All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.*

Jan 2021