

Luer Cap

IMM Luer Caps design for male or female luer connectors has different options for user preference:

- Vented or Non-vented
- Color options
- EtO compatible unless otherwise noted



Luer Cap						
P/N	Description	Use on	Vent	Color	Material	Note
A0003	Cap for female luer lock, non-vented	FLL	Non-vented	Blue	HDPE	
A0004	Cap for female luer lock, non-vented	FLL	Non-vented	Original	HDPE	
A0096	Cap for female luer lock, vented hole	FLL	Vented	Green	HDPE	
A0103	Cap for female luer lock, non-vented	FLL	Non-vented	White	HDPE	
A0123	Cap for female luer lock, non-vented	FLL	Non-vented	Transparent	PC	
A0124	Cap for female luer lock, non-vented	FLL	Non-vented	Green	HDPE	
A0125	Cap for female luer lock, non-vented	FLL	Non-vented	Red	HDPE	
A0130	Cap for female luer lock, vented	FLL	Vented	White	HDPE	
A0190	Cap for female luer lock, vented	FLL	Vented	Original	LDPE	
A0208	Cap for female luer lock, non-vented	FLL	Non-vented	Yellow	HDPE	
A0223	Cap for female luer lock, vented	FLL	Vented	Original	LDPE	
A0283	Cap for female luer lock, vented	FLL	Vented	blue	HDPE	
A0362	Cap for female luer lock, vented hole	FLL	Vented	White	HDPE	
A0425	Cap for female luer lock, vented, cover	FLL	Vented	Original	PP	
A0614	Cap for female luer lock, vented	FLL	Vented	White	HDPE	
A0615	Cap for female luer lock, vented	FLL	Vented	blue	HDPE	

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.

Jan 2021

Luer Cap						
P/N	Description	Use on	Vent	Color	Material	Note
A0365	Male/Female adapter, cap	FLL/MLL	Vented	Original	HDPE	
A0569	Male/Female adapter, cap	FLL/MLL	Vented	Yellow	HDPE	
A0601	Male/Female adapter, non-vented, cap	FLL/MLL	Non-vented	Red	HDPE	
A0020	Cap for male luer slip, vented	MLL	Vented	White	LDPE	
A0021	Cap for male luer slip, vented	MLL	Vented	White	LDPE	
A0085	Cap for male luer slip, vented	MLL	Vented	blue	HDPE	
A0091	Cap for male luer slip, vented	MLL	Vented	blue	HDPE	
A0101	Cap for male luer slip, vented	MLL	Vented	blue	HDPE	
A0133	Cap for male luer slip, vented	MLL	Vented	Original	PP	
A0194	Cap for male luer slip, vented	MLL	Vented	Original	LDPE	
A0214	Cap for male luer slip, vented	MLL	Vented	Orange	LDPE	
A0235	Cap for male luer slip, vented	MLL	Vented	Original	PP	
A0236	Cap for male luer slip, vented	MLL	Vented	Yellow	LDPE	
A0351	Cap for male luer slip, vented	MLL	Vented	Red	LDPE	
A0352	Cap for male luer slip, vented	MLL	Vented	Original	PP	
A0576	Cap for male luer slip, vented	MLL	Vented	White	LDPE	
A0577	Cap for male luer slip, vented	MLL	Vented	Original	LDPE	

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.

Jan 2021