



Tube Clamp

IMM Tube Clamp has different options for selection

- Fit wide range OD
- Different operation option
- Various color option
- EtO compatible unless otherwise noted

Subgroup	P/N	Description	Fit OD (mm)	Color	Material
Slide	A0057	Slide clamp, small, open	2.4	White	HDPE
Slide	A0090	Slide clamp, small, open	2.4	Blue	HDPE
Slide	A0131	Slide clamp, small, open	2.4	White	ABS
Slide	A0176	Slide clamp, small, open	2.4	Yellow	HDPE
Slide	A0334	Slide clamp, small, open	2.4	Blue	ABS
Slide	A0353	Slide clamp, small, open	2.4	Purple	HDPE
Slide	A0575	Slide clamp, small, open	2.4	Red	ABS
Slide	A0012	Slide clamp, small, closed	2.4	White	HDPE
Slide	A0013	Slide clamp, small, closed	2.4	Blue	HDPE
Slide	A0014	Slide clamp, small, closed	2.4	Green	HDPE
Slide	A0015	Slide clamp, small, closed	2.4	Yellow	HDPE
Slide	A0191	Slide clamp, small, closed	2.4	Red	HDPE
Slide	A0327	Slide clamp, small, closed	2.4	Green	HDPE
Slide	A0316	Slide clamp, small, closed	2.4	Purple	HDPE
Slide	A0132	Slide clamp, small, closed	2.4	White	ABS
Slide	A0224	Slide clamp, small, closed	2.4	Green	ABS
Slide	A0568	Slide clamp, small, closed	2.4	Red	ABS
Slide	A0089	Slide clamp, closed	4.1	Blue	HDPE
Slide	A0110	Slide clamp, closed	4.1	White	HDPE
Slide	A0225	Slide clamp, closed	4.1	Green	HDPE

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.

April 2023

Subgroup	P/N	Description	Fit OD (mm)	Color	Material
Slide	A0334	Slide clamp, closed	4.1	Purple	HDPE
Slide	A0449	Slide clamp, big, closed	4.1	Blue	PP
Pinch	A0289	Pinch clamp, small	4.1	Purple	PP
Pinch	A0370	Pinch clamp, small	4.1	Green	PP
Pinch	A0374	Pinch clamp, small	4.1	White	PP
Pinch	A0375	Pinch clamp, small	4.1	Blue	PP
Pinch	A0416	Pinch clamp, small	4.1	Cream	PP
Pinch	A0417	Pinch clamp, small	4.1	Black	PP
Pinch	A0307	Pinch clamp, mini	2.4	White	PP
Pinch	A0464	Pinch clamp, mini	2.4	Green	PP
Pinch	A0602	Pinch clamp, mini	2.4	Blue	PP
Pinch	A0629	Pinch clamp, small	2.4	Original	POM
Pinch	A0529	Pinch clamp, big, adjustable	9.8	White	PP
Pinch	A0558	Pinch clamp, big, adjustable	9.8	Blue	PP
Roller	S10323	Roller clamp, OD4.1mm	4.1	Original	ABS

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.

April 2023