

## Male Luer Connector & Spin Nut

IMM Male Luer Connectors are available in a wide range of configurations to meet different market requirements:

- Lock type or Slip type
- For different tube OD
- With/Without Wing
- Different spin nut for user
- EtO compatible unless otherwise noted



### Male Luer Connector and Spin Nut

P/N	Description	Tube OD	Type	Wing	Color	Material	Note
A0041	Male luer lock ID2.0mm	2.0	Lock	No	Transparent	PVC	
A0412	Male luer slip ID2.0mm	2.0	Slip	No	Transparent	PVC	D type
A0076	Male luer lock ID2.2mm	2.2	Lock	No	Transparent	PVC	
A0032	Male luer lock ID2.4mm	2.4	Lock	No	Transparent	PVC	
A0129-2B	Male luer lock ID2.4mm	2.4	Lock	Wing	Transparent blue tint	PVC	
A0175	Male luer lock ID2.4mm	2.4	Lock	Wing	Transparent	ABS	
A0238	Male luer lock ID2.4mm	2.4	Lock	No	Transparent	PVC	
A0247	Male luer lock ID2.4mm	2.4	Lock	Wing	Transparent	ABS	
A0367	Male luer lock ID2.4mm	2.4	Lock	No	Transparent brown tint	PVC	
A0528	Male luer lock ID2.4mm	2.4	Lock	Wing	Transparent	PVC	
A0059	Male luer slip ID2.4mm	2.4	Slip	No	Transparent	PVC	A type
A0173	Male luer slip ID2.4mm	2.4	Slip	No	Transparent	PVC	B type

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@imm.com.tw

*All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.*

### Male Luer Connector and Spin Nut

P/N	Description	Tube OD	Type	Wing	Color	Material	Note
A0271	Male luer slip ID2.4mm	2.4	Slip	No	Transparent	PVC	D type
A0466	Male luer slip ID2.4mm	2.4	Slip	No	Transparent	PVC	D type
A0598	Male luer lock ID2.6mm	2.6	Lock	Wing	Transparent lavender tint	PC	ETO/Gamma
A0034	Male luer lock ID2.8mm	2.8	Lock	No	Transparent	PVC	
A0138-2B	Male luer lock ID2.8mm	2.8	Lock	Wing	Transparent blue tint	PVC	
A0436	Male luer lock ID2.8mm	2.8	Lock	No	Transparent	PVC	
A0288	Male luer slip ID2.8mm	2.8	Slip	No	Transparent	PVC	A type
A0422	Male luer slip ID2.8mm	2.8	Slip	No	Transparent	PVC	B type
A0348	Male luer lock ID3.0mm	3.0	Lock	No	Transparent	PVC	
A0304	Male luer slip ID3.0mm	3.0	Slip	No	Transparent	PVC	D type
A0277	Male luer lock ID3.2mm	3.2	Lock	Wing	Transparent	PVC	
A0521	Male luer lock ID3.2mm	3.2	Lock	Wing	Transparent	PVC	
A0399	Male luer slip ID3.2mm	3.2	Slip	No	Transparent	PVC	B type
A0424	Male luer slip ID3.2mm	3.2	Slip	No	Transparent	PVC	D type
A0028	Male luer lock ID3.25mm	3.25	Lock	Wing	Transparent	ABS	
A0443	Male luer slip ID3.3mm	3.3	Slip	No	Transparent	PVC	B type
A0409	Male luer slip ID3.3mm (2.4mm)	3.3	Slip	No	Transparent	PVC	B type
A0035	Male luer lock ID3.5mm	3.5	Lock	No	Transparent	PVC	
A0137-2B	Male luer lock ID3.5mm	3.5	Lock	Wing	Transparent blue tint	PVC	
A0205	Male luer lock ID3.5mm	3.5	Lock	No	Transparent	PVC	
A0233	Male luer lock ID3.5mm	3.5	Lock	Wing	Transparent lavender tint	PC	ETO/Gamma
A0061	Male luer slip ID3.5mm	3.5	Slip	No	Transparent	PVC	A type
A0179	Male luer slip ID3.5mm	3.5	Slip	No	Transparent	ABS	C type
A0195	Male luer slip ID3.5mm	3.5	Slip	No	Transparent	PVC	B type
A0217	Male luer slip ID3.6mm	3.6	Slip	No	Transparent	PVC	B type
A0412	Male luer slip ID3.8mm	3.8	Slip	No	Transparent	PVC	B type
A0033	Male luer lock ID4.1mm	4.1	Lock	No	Transparent	PVC	

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: [service@imm.com.tw](mailto:service@imm.com.tw)

*All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.*

Jan 2021

### Male Luer Connector and Spin Nut

P/N	Description	Tube OD	Type	Wing	Color	Material	Note
A0136-2B	Male luer lock ID4.1mm	4.1	Lock	Wing	Transparent blue tint	PVC	
A0234	Male luer lock ID4.1mm	4.1	Lock	Wing	Transparent lavender tint	PC	ETO/Gamma
A0582	Male luer lock ID4.1mm	4.1	Lock	No	Transparent	PVC	
A0060	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	PVC	A type
A0178	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	ABS	C type
A0196	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	PVC	B type
A0272	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	PVC	D type
A0434	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	PVC	D type
A0451	Male luer slip ID4.1mm	4.1	Slip	No	Transparent	PVC	C type
A0621	Male luer slip ID4.1mm (small outlet)	4.1	Slip	No	Transparent	PVC	B type
A0002	Male luer lock ID4.4mm	4.4	Lock	No	Transparent	ABS	
A0426	Male luer lock ID4.4mm	4.4	Lock	No	Transparent	PVC	
A0275	Male luer lock ID4.8mm	4.8	Lock	No	Transparent	PVC	
A0536	Male luer lock ID5.3mm	5.3	Lock	No	Transparent	PVC	
A0218	FLL to MLL	na	Lock	No	Transparent lavender tint	PC	ETO/Gamma
A0543	FLS to MLS	na	Slip	No	Transparent	PVC	D type
A0535	MLL to MLS	na	Slip	No	Transparent	PC	
A0005	Spin nut, A type	na	Lock	No	Transparent	ABS	
A0115	Spin nut, A type	na	Lock	No	Transparent	PC	
A0209	Spin nut, A type	na	Lock	No	White	ABS	
A0287	Spin nut, A type	na	Lock	No	Transparent lavender tint	PC	ETO/Gamma
A0172	Spin nut, B type	na	Lock	No	Transparent	ABS	
A0423	Spin nut, B type	na	Lock	No	Transparent lavender tint	PC	ETO/Gamma
A0177	Spin nut, C type	na	Lock	No	Transparent	ABS	
A0450	Spin nut, C type	na	Lock	No	Transparent lavender tint	PC	ETO/Gamma

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: [service@imm.com.tw](mailto:service@imm.com.tw)

*All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.*

Jan 2021