

Multiway Connectors

IMM provide multiway connectors for different application design:

- Two (2), three(3), four(4) way to one (1)
- Tube or luer distal end
- Different tube OD

Material EtO compatible unless otherwise noted



P/N	Description	Way	Tube OD	End	Color	Material	Note
A0058	Mini-Bi-RMLS OD2.4mm	2	2.4	Slip	Transparent	PVC	
A0075	Mini-Bi-RMLS OD2.2/2.4mm	2	2.4	Slip	Transparent	PVC	
A0114	Mini-Bi-RMLS OD2.4mm-long	2	2.4	Slip	Transparent	ABS	
A0152	Mini-Bi-RMLS OD2.4mm	2	2.4	Slip	Transparent	ABS	
A0232	Mini-Bi-RMLS OD2.2/2.4mm-long	2	2.4	Slip	Transparent	ABS	
A0414	Mini-Bi-RMLS OD2.4mm-long	2	2.4	Slip	Transparent	PVC	
A0383	Bi-way connector OD3.2mm	2	3.2	Tube	Transparent	PVC	
A0025	Bi-way connector OD3.5mm	2	3.5	Tube	Transparent	PVC	
A0026	Bi-way connector OD4.1mm	2	4.1	Tube	Transparent	PVC	
A0581	Bi-way connector OD6.4mm	2	6.4	Tube	Transparent	PVC	
A0364	Bi-way connector OD6.7mm	2	6.7	Tube	Transparent	PVC	
A0027	Tri-way connector OD4.1mm	3	4.1	Tube	Transparent	PVC	
A0427	Qua-way RMLS OD2.4mm	4	2.4	Slip	Transparent	PVC	

Innovative Medical Manufacturing Company

107, Lane 181, Sec. 1., Yongzhen Rd., Zhunan Township, Miaoli 35057, Taiwan R.O.C.

TEL: +886-37-620236 FAX: +886-37-620239 E-mail: service@immc.com.tw

All IMM components are supplied bulk, non-sterile for further processing before clinical use. The buyer is responsible for the appropriateness of the final application.

Jan 2021